

# Hydrocolloid Dispenser by TCS



The Ideal machine for small to medium size labs.

## 2 & 3 Gallon

- Easy to use and maintain
- Double Boil Out Tank
- Automatic melting & cooling
- No need for water connections

*Available by Special Order!*

**10 Gallon Unit**

High Production Labs



Made in The USA  
5 year warranty

### Built Tough! Machines built in 1966 still in use!

One of the only double boiler tanks on the market, and **THE ONLY DOUBLE BOILER THAT USES ANTI-FREEZE** instead of circulating water. This engineering eliminates the need to connect the machine to a water source and to drain water periodically. The Anti-Freeze never needs to be changed, simply check the level of anti-freeze every few months and maintain at specified level.

The machine pays for itself by doubling the life of your hydrocolloid material. The double boiler does not allow the material to burn, a significant saving for your laboratory.

The double boiler spares the hydrocolloid paste from sticking to the sides of the tank. This means that the mixer turns periodically instead of constantly (like other duplicating machines), saving on electricity and needless use of mechanical parts.

### 2 Gallon

**Specifications:**

22 inches high  
Tank diameter 8"  
Width 13" (including control box / pour spout)  
Approximately 20 lbs.  
Available in 110V or 220V

### 3 Gallon

**Specifications:**

21 inches high  
Tank diameter 12"  
Width 18" (including control box / pour spout)  
Approximately 25 lbs.  
Available in 110V or 220V

### 10 Gallon

**Specifications:**

33 inches high  
Tank diameter 14.5"  
Width 19" (including control box / pour spout)  
Approximately 75 lbs.  
Available in 110V or 220V

Distributed by: